

BOSTON BAND

Manufactured by Boston Orthotics & Prosthetics, Avon, MA 02322

Indications for use: The Boston Band is intended to passively hold prominent cranial regions of an infant's skull to improve cranial proportion and symmetry in infants from three to eighteen months of age, with non-synostotic positional plagiocephaly, including infants with plagiocephaly, brachycephaly, and scaphocephaly. The device is also indicated for adjunctive use for infants from three to eighteen months of age whose synostosis has been surgically corrected, but who still have moderate to severe cranial deformities.

Contraindications: Untreated craniosynostosis or hydrocephalus.

Precautions: This device may cause skin irritation or breakdown. Contact your treating orthotist in the event of an adverse event.

Instructions for use: Wear & Care Guide to be provided to caregiver.

Caution: Federal law restricts this device to sale by or on the order of a physician and treating orthotist.

Revised: 2022-12-1



BOSTON BAND

For a digital copy of the Wear & Care Guide and other helpful information on the Boston Band, scan the code below.



